SNHELLVERASCHER

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RAPID INCINERATOR

up to 999°C

for single and series incineration of solids
and liquid substances

Rapid incinerators are used for time-saving single or series incineration of solid and liquid substances. The aim is often the analysis of substances and liquids for their organic and inorganic components. Due to the two-story structure, preheating, drying, fuming, and finally, ashing can co-occur in only one unit. The operating temperature of 999°C is reached in a few minutes.

The upper platform has 8 recesses with a diameter of 35 mm which can accommodate porcelain crucibles of up to 50 mm diameter for drying, pre-heating and fuming-off. The uninhibited exposure to oxygen speeds up the incineration process. The inclined position of the crucibles in the lower incinerating bed (which has a length of 400 mm) allows excellent observation and access to the contents during operation.

Equipment for SVII:

The upper part and the main housing are made of high-grade stainless-steel; the housing comes with an orange finish (RAL 2002). The device has a built-in, 120-minute timer with an audible signal for process end and permanently ON position. The appliance is fitted with a safety switch which cuts off the mains supply to the heater immediately when the service flap is opened. The unit is delivered ready for connection with approx. 1.8 m power cable and Schuko plug.

Equipment for SVR/E:

As SV/II but with an additional electronic temperature regulation by means of a performance control.

Equipment for SVD95:

As SV/II, but with an additional PID temperature regulator up to 999°C, this is controlled by a built-in, NiCr-Ni thermo couple. The pre-selected temperature is shown on the digital display with a constant ACTUAL temperature.

Best-Nr.	SVII	SVR/E	SVD95	Order No.
Spannung	230 - 240 V, 50-60 Hz			voltage
Leistung	2500 Watt			performance
Breite x Tiefe	450 mm x 180 mm			width and depth
Höhe	310 mm			height
max. Tiegel-Ø	50 mm			max. crucible diameter
max. Tiegelhöhe	50 mm			max. crucible height
Gewicht	7,4 kg	8,0 kg	7,8 kg	weight
Best.-Nr. Ersatzheizelement	EH 400	EH 400	EH 95	Order No. spare heating element

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The right to make modifications serving further technical development, is reserved.

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Further Fabrications

High temperature TITANIUM hot plate up to 600°C - CERAN heating plates to 500°C - Precision hotplates - Precision hotplate made of cast iron up to 450°C - Stand hotplates with supporting rod - Multipurpose, extremely deep, heating bath - Sandbaths - Series heating units - Multipurpose heating baths - Water- and Oil baths - Development - Special orders.